

# Dektak-xp Stylus Profilometer Operating Procedure

Effective Date: 09/03/2012

Author(s): Jiong Hua

Phone: 402-472-3773

Email: [jhua2@unl.edu](mailto:jhua2@unl.edu)



# 1 Introduction

## 1.1 Key Words

Stylus profilometer

## 1.2 Purpose

This document provides instructions for operating the Dektak-xp stylus profilometer system, which provides high resolution thin film profiling capability. Use of this tool requires the understanding of the fundamentals of stylus type profilometer knowledge.

## 1.3 Applicability

### 1.3.1 Locations

The tool is located at **clean room of nanoscience research center.**

### 1.3.2 Safety

The safety concern is the system uses high voltage electricity and high pressure gases.

### 1.3.3 Restrictions and Limitations

- **Must be a qualified user of clean room.**
- **The limitation for the substrate is that it must not damage the system. Sample heights must NOT exceed 50mm**



## 1.4 Restrictions on Working Alone



- Normal working hours are from 8am to 6pm M-F.
- Working alone is permitted with completion of an orientation to this written procedure and hands-on training from the specialist.
- Assistance from the specialist is available during working hours only. If an error occurs during off-hours, record the error in the Logfile and send an email to the specialist. **Do not try to fix or adjust anything by yourself.** Tool will be checked in the following work day. User will be notified when sample left in chamber is available for pickup.
- Problems with equipment malfunctions, breakage, etc. should be reported to the specialist and recorded in the tool Logfile. **Again do not try to fix or adjust anything by yourself.**
- In case of emergency, press the “EMO” button on the left side of the system. Contact the specialist right away.
- For any emergency involving injuries, fire, chemical spills, etc., call **911.**



# 2 Preparations




- Receive this procedure from the specialist.

### 3 Execution (Step-by-step work breakdown)

Step #	Alerts	Action
1.		<p><b>Checking status of the system before starting:</b></p> <ul style="list-style-type: none"> <li>• Check and ensure the compressed air gauge is at ~10 psi.</li> <li>• Check and ensure the main power is off.</li> <li>• Check and ensure that the computer is running.</li> </ul>
2.		<p><b>Starting the system:</b></p> <ul style="list-style-type: none"> <li>• Turn on the main power of the system by pressing “I” on power control unit.</li> <li>• Logon to the computer.</li> <li>• Open and fill out the logfile on the desktop.</li> <li>• Launch the “Vision64” software on desktop. The system will initialize the stage when the software starts. You need to confirm the initialization in both XY direction and Theta.</li> <li>• Wait until the stage stopped and user interface appears on the screen. To get the best result, allow the system to warm up about <b>15 minutes</b> before take a measurement.</li> </ul>
3.		<p><b>Loading the sample:</b></p> <p><b>Note:</b> Sample heights must NOT exceed <b>50mm</b>.</p> <ul style="list-style-type: none"> <li>• Open the door of the environment chamber.</li> <li>• Click “Unload sample” button  on the tool bar above the live video display window to move the stage to loading sample position.</li> <li>• Wait until the stage fully stopped, put sample at the center of the stage.</li> <li>• Turn on vacuum switch at right corner of the floating table to fasten the sample on the stage.</li> <li>• Click “Load sample” button  on the tool bar above the live video display window to move the stage back to its home position.</li> <li>• Use XY control panel and Theta control panel in the instrument control window to move the area of interest on the sample roughly beneath the stylus.</li> </ul>

<p>4.</p>	<p><b>Locating the area of interest:</b></p> <ul style="list-style-type: none"> <li>• Click “Tower down” button  on the tool bar above the live video display window to lower the tower assembly down to bring stylus to contact level.</li> <li>• Choose suitable intensity (normally <b>5%</b>) for CCD camera at the low right corner of the main window. The live video display window now shows the detail of the sample surface.</li> <li>• Use XY control panel and Theta control panel in the instrument control window to move the area of interest on the sample right beneath the stylus. <b>Make sure the stylus won't touch the sample surface during this procedure.</b></li> <li>• <b>Right click on live video display window and click on “Lower stylus” in the pop-up menu.</b></li> </ul> <p><b>Note:</b> The current position of the stylus is the starting point of the scanning which will act along Y axis from bottom to top.</p>
<p>5.</p>	<p><b>Measurement setup:</b></p> <p>In measurement setup window, choose suitable parameters:</p> <ul style="list-style-type: none"> <li>• Scan type: “Standard scan” (for 2-D scanning).</li> <li>• Range: Choose one most close to yet bigger than the maximum height of your sample. If unknown, start from a bigger one (524um).</li> <li>• Scan profile: Choose one fits your sample profile. If unknown, choose “Hills and valleys ”.</li> <li>• Stylus type: Choose “radius 50nm” unless other acclaimed.</li> <li>• Stylus force: Enter between 1mg to 15 mg. <b>Use 10mg for most samples.</b></li> <li>• Length: Scanning distance (from 50um to 55,000um)</li> <li>• Duration: Scanning time. Adjust it so that the Resolution is no less than the tip size (<b>50nm</b> unless other acclaimed).</li> <li>• Sampling points and speed are calculated by the system itself based on “Length” and “Duration”. Users can not change them directly.</li> </ul>

6.		<p><b>Taking measurement:</b></p> <ul style="list-style-type: none"> <li>• Close the door of the environment chamber.</li> <li>• Press “Measurement” button  on the tool bar of the main window to start scanning. The window will automatically switch to Data Acquisition view.</li> <li>• Scan can be stopped at any time by pressing “cancel” button  on the tool bar of the main window. You may change the measurement set such as distance, range, etc.</li> <li>• If acquired data reaches the maximum of the range, stop scan and change the range to bigger value.</li> </ul>
7.		<p><b>Data Analysis:</b></p> <ul style="list-style-type: none"> <li>• Wait until the scan finishes. The window will automatically switch to data analysis view.</li> <li>• Applying software leveling to set the reference at zero: <ol style="list-style-type: none"> <li>1. Click “Terms Removal (F-Operator)” in the analysis tool box to add it to data analyzer tree.</li> <li>2. Position “R” (reference) and “M” (Measurement) cursors as far apart as possible along the base line.</li> <li>3. Right click on data analysis display and select “two point linear fit” from “Data leveling settings” pop-up menu. Then click “calculate” button.</li> <li>4. The profile trace re-plots and levels with the “R” and “M” cursor intercepts at zero.</li> </ol> </li> <li>• Step analysis: <ol style="list-style-type: none"> <li>1. Move “R” to reference zero and “M” to the step area.</li> <li>2. Drag and expend the width of “R” and “M” cursor.</li> <li>3. The result appears in a box at the bottom of data analysis window.</li> </ol> </li> </ul>

8.		<p><b>Unloading the Sample:</b></p> <ul style="list-style-type: none"> <li>• Press “Measurement Setup” button in the tool bar of the main window to bring the view back to measurement view.</li> <li>• Press “Tower Home” button  on the tool bar above the live video display window to bring the tower assembly all the way up to its home position.</li> <li>• Click “Unload sample” button  on the tool bar above the live video display window to move the stage to loading sample position.</li> <li>• Wait until the stage fully stopped; open the door of the environment chamber.</li> <li>• Turn off vacuum switch at right corner of the floating table to release the sample on the stage and take the sample out.</li> <li>• Click “Load sample” button  on the tool bar above the live video display window to move the stage back to its home position.</li> <li>• Close the door of the environment chamber.</li> </ul>
9.		<p><b>Switch off Procedures:</b></p> <ul style="list-style-type: none"> <li>• Close the “Vision64” software.</li> <li>• Fill out the log file with the corresponding parameters. <b>DO NOT log off the computer.</b></li> <li>• Switch off the main power of the system by pressing “0” on power control unit.</li> </ul>
10.		<p><b>Clean up all samples, pens, and notebooks from the area.</b></p>
End		<p><b>End of Procedure</b></p>

## 4 Post-Performance

### 4.1 Recordkeeping

Completely fill out the logfile.

### 4.2 Feedback

Report any unusual or problematic behavior of the setup by contacting the specialist.

## **5 References**

### ***5.1 Technical References***

- **Dektak-xp** User Manual

## **6 User Access Level**

Normal User – Requires specialist to be present

Expert User – Does not require specialist to be present